

Go from cells to PCR in less than 10 mins!

Nippon Genetics Europe is proud to introduce a novel reagent for the lysis of cells and the release of genomic DNA. No spin columns or extractions, just a 10 minute incubation on a thermal cycler and the DNA is ready for the PCR process.

DNAreleasy has been successfully used to lyse the following:

Human Buffy coat, Human fibroblasts tissue culture cells, Human sperm

Mouse & rat tissue culture cells, Mouse tail.

Various plants (leaf, flower and flour): Arabidopsis, Cabbage, Oilseed Rape, Oats, Soya & Maize.

Flea, Drosophila and others like yeast, Gram +ve bacteria, Gram -ve bacteria, Viruses and many more ...

Protocol:

1) Mix cells with 20 μ l of DNAreleasy. It is important that the lysis solution covers the cells completely.

2) Overlay with mineral oil if necessary (for Thermal Cycler without heated lid)

3) Place in a Thermal Cycler

Step 1: 75°C for 5 mins

Step 2: 96°C for 2 mins

Step 3: 20°C hold

After the lysis, a part or all of the DNAreleasy-lysate can be used directly in PCR. DNAreleasy can make up to 40% of most PCR mixes. Or it can be stored at -20°C for future use.




Ordering information


LS01	DNAreleasy; 10 preparations
LS02	DNAreleasy; 50 preparations

DISTRIBUTED BY







Applications: Mouse Tails and Ears




Competitor V
Overnight incubation



Pro K digestion
Overnight incubation



DNAreleasy (7 min.)




Homemade method

Data kindly provided by Dr. John Foskett, Kapa Biosystems, Cape Town; South Africa.

Comparison of 4 different lysis methods. The lysate was used in the amplification of 8 different amplicons by PCR using KAPA2G Robust Hot Start ReadyMix (Kapa Biosystems, Cape Town, South Africa).

Applications: HT-29 colon carcinoma cells




1	DNA marker
2	DNAreleasy lysate, cell culture of HT-29 colon carcinoma cells, amplified with Kapa 2G Robust Hot Start
3	Plasmid DNA, amplified with Kapa 2G Robust Hot Start
4	Plasmid DNA and cell culture lysate, amplified by Kapa 2G Robust Hot Start
5	Negative control
6	DNAreleasy lysate, cell culture of HT-29 colon carcinoma cells, amplified with Kapa 2G Robust Hot Start (a different primer set was used)
7	Plasmid DNA, amplified with Kapa 2G Robust Hot Start (a different primer set was used)
8	Plasmid DNA and cell culture lysate, amplified by Kapa 2G Robust Hot Start (a different primer set was used)
9	Negative control


Data kindly provided by Ms. Gudrun Wahlström, Faculty of Medicine, University of Turku, Finland.

Statement: *"This kit worked very well without any optimization. It is a very simple procedure and the combination with the Kapa 2G Robust Hot Start enzyme leads to superior results"*.

Applications: Arabidopsis and oats: Seeds, leaves and buds



1. Arabidopsis (actin primers) 445 bp fragment
Lane 1 bud - 2 bud - 3 bud - 4 bud - 5 leaf - 6 leaf - 7 leaf - 8 leaf - 9 +ve control - 10 -ve control - 11 size ladder



2. Oats (actin primers) 380 bp fragment
Lane 1 -ve control - 2 seed - 3 seed - 4 seed - 5 leaf - 6 leaf - 7 leaf - 8 leaf - 9 +ve control - 10 size ladder